From: Steve Ross <u>Skross@anl.gov</u> Phone 630 252-9510 fax 630 252 9350 401-c1252 APS/ASD/Electrical Systems Group

### **ABSTRACT:**

Description of electrometer amplifiers applied to PIN diode based x-ray beam position and intensity monitors. This document is work-inprogress, suggestions welcome.

APS/ASD/ESG has developed an electrometer amplifier sensitive in the range of picoamps to microamps. This low cost amplifier has a high (10<sup>6</sup>) dynamic range, and also a bandwidth of up to 1600 samples/sec. Its noise floor is in the range of 100 pA RMS (1600 Hz). In one 9" x 5" package, called an ADCMOD2 module, there are 4 such electrometers. The main application at the APS is to monitor the output of 4 detectors, commonly PIN diode detectors, for x-ray beam position monitoring. Data from these modules travels over fiber optic links to a VME card set. The data then can be read from the VME bus. We believe that a basic system would read 16 PIN diode current sources, but there are several variations. We work with the beamline personnel to get a system up and running.

PIN diode Users Guide 12-17-02.doc Rev date: 12/14/2002

**CONTENTS:** This memo describes electronics developed to read currents in the range of picoamperes to microamperes. The primary application is to serve as an electrometer for various arrangements of x-ray detectors. For example this electronics is being used to measure the current from PIN photodiodes, which in turn serve as x-ray beam position and intensity monitors. (X-BPM).

This memo is organized as follows:

(1) We give background information about why this electronics was developed and where it was initially applied.

- (2) We give a description of the x-ray detector typically used a four quadrant array of PIN diodes.
- (3) We give a specification for the electrometer, based on design and on data taken at sector 32.
- (4) We describe the data acquisition system based on a fiber optic link back to a VME card set. We show what VME cards are necessary to make the system work and discuss some of the options. A basic system is defined to be one running 16 PIN diodes.
- (5) We discuss software requirements. The electronics can be controlled remotely, from the VME bus. One-time output setup scripts are described. We describe how to read the returned data.
- (6) We list a number of issues that typically must be addressed to get the system at a particular beamline at the APS. We work with you to install a system and we provide almost all parts. Costs are discussed.
- (7) We give information useful to running the system. This section is continually being edited, to be clearer, and as we gain more experience with users running registered mode.
- (8) Schematic information. This section will grow into a description of where all the design files are documented. It will mainly be of use to ESG personnel.
- (9) Future features. This discussion will be never ending.
- (10) APPENDIX. For now this is where I add information mainly relevant to the digital storage scope mode, and sector 32 operation. It will however help me to go over this mode with beamline users who may want to move over to it.

(1) **BACKGROUND.** This electronics was originally developed per the specifications of ComCAT sector 32 at the APS [REF: Steve Wasserman, Kevin D'Amico]. This work in turn was inspired by an RSI paper discussing results from an earlier electrometer design. A copy of this electronics was installed at COMCAT in early 2000, with the primary application of EXAFS measurements. Their electronics collects 12 PIN diode signals, plus the signal from the encoder on the x-ray monochromater. As the energy is swept, plotting one versus the other gives an EXAFs curve. Thus this system works in "digital storage scope mode", as contrasted with "register mode". (These modes of operation will be further discussed below.)

(2) PIN DIODE DETECTOR. The detector arrangement was based on designs at sector 19 [REF: R.W.Alkire, G. Rossenbaum, G.Evans J. Sync Radiation (2000) 7 p.61]). Randy gave me a copy of this article in pdf format. The basic idea is that the monoenergetic x-ray beam passes through a hole in the circuit board holding the PIN diodes, and then on through a thin metal foil, and then on to the experiment. The backscattered fluorescence from the metal foil is detected by the four diodes arranged as a quadrant. As the beam moves closer to one diode, its signal strength will increase. Standard centroiding or difference over sum calculations can then show beam position. I believe that only Randy at sector 19 has actually moved the quad PIN detector in a calibrated way, and this was done without using this set of electronics. But based on his

work, and the specifications of this system, it is likely that the position of the beam can be estimated with sub-micron precision.

The red wire of the PIN diode connects to the signal input of the amplifiers; the black wire connects to the ground. The PIN diodes have a common (analog) ground at the printed circuit board. Four pairs of wires go from the amplifiers to the four PIN's.

Randy works with companies which supply the foil xxxxname?. One issue is that the foil needs to be without "wrinkles", that it be flat. Otherwise, as the APS SR x-ray beam moves across the foil, the signal will change simply from the ups-and-downs of foil itself. The vendor will stretch the foil over a beveled washer, to make it tight like the head of a drum. With micron thick foils, beamlines report current levels of 100's of nanoamps, a good match to the electrometers described here.

Randy also makes the following comments about use of this type detector [REF: email 12/11/02]:

1. All metal foils are pressure sensitive but usually compare well to the response one gets from an ion chamber i.e., response versus energy. For kapton backed foils, the response may differ from an ion chamber due to the multiple elements involved in both fluorescence and scattering.

2. Separating the diodes from the foil is a good idea but not essential to performance. It does allow one to turn the BPM on and off by simply removing the foil. Also if the diodes are moved and not the foil during calibrations the uniformity of the foil becomes less of an issue.

3. These diodes are light sensitive and that means ANY light, including light from ion gauges if they are nearby.

4. Foils can become contaminated with carbon deposits over time and may require replacement. Not an issue for operations of less than a year.

5. This device is only for monochromatic operations, not white light.

The mechanical design for the PIN diode holder has since undergone several design changes. These were done for the convenience of other APS CAT's. Some beamlines have unique mechanical constraints. **The point is that in short order we can work with the beamline users to develop custom printed circuit boards for this mounting.** Three examples of such printed circuit boards are shown. [FIG 1 A,B,C,D,E]



FIG 1A This is the basic printed circuit board type design (DIODE5) holding the four UDT 1 cm2 PIN diodes. It is about a small as possible with these diodes. The center hole in milled out to about 0.5 inch diameter. This circuit board fits into a KF-40 or KF-50 size vacuum system.





FIG 1 B,C A photograph, and a drawing of the diode6 design. This printed circuit board fits into a KF50 cross and positions the diodes in the center. Standoffs hold a second circuit board approximately 1 cm. above it and this circuit board holds the foil. This is a very simple one-piece design for medium vacuum conditions.



FIG 1D Photograph of DIODE1 layout. Again this is just another variation of how the diodes can be arranged. Red wires signal, black wire are grounded together on the PCB.





FIG 1E this is the design for the diode detectors at sector 32. (DIODE1) The diodes are on a circuit card, which can slide into a cardholder type mount. At this sector, the metal foils can be changed remotely, so that different foils can be inserted for different EXAFs scans

## (3) SPECIFICATION FOR THE ELECTRONICS.

We have two electrometer modules based on an integrated charge amplifier and based on a voltage amplifier. Almost all the discussion here relates to the charge amplifier. It is the version currently going out onto the APS storage ring beamlines. We anticipate the use of the voltage amplifier when the signal currents exceed about 5 uA. Grounding for optimum noise performance is an open issue to be addresses at each beamline.

**CHARGE AMPLIFIER**. The analog characteristics of the electrometer are based on the Texas Instruments Inc. (formerly Burr Brown Inc.) ddc112 integrated circuit. To see its specifications, go to www.ti.com and search on the ddc112 chip. A photograph of the printed circuit board is given. (FIG 2) All the modes and options discussed in the data sheet are implemented in our system. However because ultimately the integrated circuit uses a 20 bit analog-to-digital converter (ADC) we find that most options, most "bells and whistles" are not needed. Twenty bits is a dynamic range of a million. When run in a somewhat standard manner (using the external 220 pf capacitor for maximum rated full well), we observed the following in digital storage scope mode:

•4 Channels - currents from 4 unbiased PIN diodes processed in each module

•Maximum readout rate 1600 samples/sec per diode [1]

•Noise of electronics alone <25 pA RMS (has not yet been optimized)

•Noise of system, including UDT S100-VL 1 cm<sup>2</sup> diode <100 pA RMS (battery), 200 pA RMS (lab supply)</li>
•Dynamic Range: 20 bit ADC, LSB = 1 pA, full = 5 uA [2]
•Typically all charge collected (one amplifier is collecting while another is processing through ADC)
•Variable gain, variable integration time, set by VME commands [3]
•Fiber optic digital readout back to/from VME
•Software controllable (see below)

Notes:

- [1] See discussion of conv command.
- [2] Will change with signal averaging.
- [3] See the range command and the conv command

We found that the twisted pair wire that connects the PIN diodes to the electrometer can be up to 2 meters in length, we have not pushed this limit yet. Thus the electrometers do not need to be located right at the detector.

Each ADCMOD2 module requires about .38 amps. The supply voltage should be in the range of 9-14 volts, I typically use +10 volts lab supply. This supply is "cleaned up" with a voltage regulator chip. The units can be powered by a +12 VDC battery, which can last for about 12 hours. (See also future-features section)

The ddc112 data sheet discusses several parameters that can be varied, the most important being gain control and its own readout mode (not to be confused with readout modes of the overall system, discussed below). Gain range control is performed by varying the feedback capacitor. (See range command.) More often used is gain control implemented by varying the integration time. Of course longer integration times slow down the overall data rate. The twoddc112 readout modes are "continuous" and "non-continuous" modes. The charge amplifier can either integrate all charge coming to it (continuous), or, if the signal level is high, it can integrate short windows, then go insensitive (discontinuous). These short windows of operation really do not speed up the sample rate however, as the internal conversions and processing still consume a fixed overhead time. A readout rate of about 1600 Hz is about as fast as things go. We mention that we have not yet pushed the limits of gain control. Eventually we will place on the circuit board feedback capacitors larger than those specified in the data sheet, with the goal of being able to measure higher current levels with acceptable linearity.

As we discussed this system with other APS users, we realized that most users actually want continuous readout of the XBPM. However about 1600 samples per second per PIN typically is too high for many beamline data systems. Data cannot be read into EPICS at this rate. We mitigate this with digital signal averaging in the PASSTH4 card.



FIG 2. Charge amplifier four channel electrometer. It is placed into a box measuring 9" x 5.5" x 2.5".

**VOLTAGE AMPLIFIER**. Another electrometer printed circuit board has been fabricated and partially tested. At present I am not spending a lot more time on this path, but can return to it in the future. The main driver to create this system was from users who wish to monitor higher levels of currents – microamps on up through milliamps. The system can work at lower currents, but when we cross back into nanoamps, the charge amplifier works very well.

Basically this amplifier consists of an electrometer grade operational amplifier (Burr Brown OPA129U). There are three gain ranges, remotely selectable. When a different gain range is selected, a micro relay switches in a different feedback resistor for the amplifier. The analog amplifiers then send signal into a 24 bit 40 KHz ADC (ADS 1252U). (FIG 3)



FIG 3. A photograph of the voltage amplifier method for electrometers. This board holds four chains of electrometer voltage amplifiers. Data is digitized, and can be packaged in a way consistent with the charge amplifiers. There are V/F converters on this board but this is not a desired development path, ADC's are better.

This hardware exists but work needs to be done to program the firmware. This hardware is mechanically compatible with the charge amplifier approach. We could build another circuit which would allow toggling between the charge and voltage amplifiers – more bells and whistles for the future.

# (4) VME BASED DATA ACQUISITION – WHAT YOU NEED FOR BASIC SYSTEM.

We discuss a basic system – what you need. We also discuss here two system readout modes, with the emphasis on the registered mode. Keep in mind that these cards really just form a VME based data acquisition system. In other words, as an option, we can input into the data stream "other" sources of data such as the timing of the APS top off signal. There is much flexibility here both in terms of hardware (connectors to input/output the data) and firmware.

#### 4.1 BASIC SYSTEM.

We can define a very basic system as supporting four PIN diodes, one ADCMOD2 card. The next level up is to have a system that supports16 PIN diodes, and this is the system we describe here. To do this the beam line needs:

- (1) Qty 4: the quad PIN diode mechanical holders. See Section 2 of this memo.
- (2) Qty 4: Four channel electrometer electronics boxes, "ADCMOD2". These will reside in the hutch, near the detectors. They are the amplifiers. (FIG 2,3)
- (3) Qty 4: terminated fiber optic cable pairs run from the hutch, to the control rack. Each ADCMOD2 board needs a pair of fibers, one transmits, one receives. Typically the fiber is bundled as 4,8,12 fibers in a cable. For example the breakout fiber (orange color) has 4 125/62.5 um fibers in it. We can provide the fiber.
- (4) Qty 1: an APS supplied VME card "FIBER2" to receive the fiber optics (FIG4) One of these cards can support four (4) ADCMOD2's.
- (5) Qty 1: an APS supplied VME card "PASSTH4", which can be read by a VME host, such as an IOC running EPICS. (FIG 5)

As discussed in section 3, the ADCMOD2 amplifiers read in the PIN diode signals, convert to digital levels, and output via a fiber.. Each FIBER2 card can support 4 ADCMOD2 cards, which in turn read 16 PIN diodes (e.g. four quad arrays). Each FIBER2 card needs the PASSTH4 card to format the data to be read over the VME backplane.



Fig 4 FIBER2 VME board



#### Fig 5. PASSTH4 VME board.

**4.2 REGISTER MODE vs. DIGITAL STORAGE SCOPE MODE**. Data can either continuously come out of the ADCMOD2 cards (registered mode), or can be output for a time, and then stop when a large memory is filled (digital storage scope mode). The emphasis in this memo is on the "register mode" of operation.

**Registered Mode**. The purpose of the register mode of operation is to take the data from the PIN diodes and get it to a location readable over the VME bus. The data comes out of all the ADCMOD2's, travels over the fiber back to the control VME rack. It is moved into a static RAM memory in the PASSTH4 card. This memory can be read over the VME bus by an A24/D16 parallel input/output (PIO) command. The data is continuously being updated. Thus the register values will continuously change as new data arrives from the electrometers.

Thus to read the signal of a given PIN diode, the VME controlling computer would access a specific VME A24 address. For example read the D16 values from hex address 0xf000, 0xf004, 0xf008, 0xf00c and you obtain the 16 bit PIN values. This would typically be a parallel input/output VME operation (PIO).

This is how the system stands now (12/2002). Two problems have now become obvious to me: (1) few beamlines can read the data OUT of the VME addresses as fast as the

electrometers can put the data IN. (The values in the registers would likely change about every 800 us.) Thus I have added signal averaging. PASSTH4 has the functionality that it sums up the signals from a given electrometer (PIN diode), and then updates the VME register at a slower rate. This averaging also improves the noise performance (at a trade off of bandwidth). The exact number of samples that are summed is TBD, not to exceeed 4096. This parameter will become more "user friendly" to change as I improve the firmware. (2) The data samples, or averaged data samples are not co-registered in time. If the VME system is not fast enough, it may read the PIN signal from diode 1 and 2 (for example) at time t1, but not get back to reading the signal from diode 3 and 4 until later. This can be fixed without much trouble but it is not done yet. I am learning user requirements as I go.

**Digital Storage Scope Mode.** The original sector 32 electronics was meant for EXAFS, and as such recorded simultaneously monochromater encoder position signals. A user-specified number of samples were recorded into a large (64MB) VME memory. (At a typical time between samples of 800 us, 1300 HZ, perhaps 50,000 samples), When the specified number of samples is collected, system stops acquiring data. The user is then expected to read all this data from the VME memory into the host computer memory.

To implement this mode of operation, we need to add two more VME cards to the system, for a total of four cards. The data from the PASSTH4 VME card discussed in section 4 is passed to a similar card, called the DC2. The purpose of the DC2 card is to handle all memory access. From the DC2, data is pumped into a commercial VSB/VME dual ported D-RAM memory. Data flows into this memory over the VSB bus. This flow does NOT interfere with VME commands on the VME backplane. Once it is set up, and told to acquire, it will handle the movement of all data to the Chrislin memory. The Chrislin memory can be read out whenever the users wishes, as it is dual ported. The intent is that the user wait until the "storage scope" type scans are complete before commencing this long transfer out.

In reality, the DC2 card appears identical to the PASSTH4 card from which it receives data, but the firmware inside the Altera chips is different. The Chrislin, Inc.VME40 memory typically costs less that \$700 for 64 MB cards. See www.chrislin.com. This company has provided such products for >12 years.

(5) SOFTWARE REQUIREMENTS. I run an entire stand-alone software system based on a PC, Microsoft Windows, and MSVC++ 6.0 code. I don't expect this to be particularly popular around the APS. I can bring it out for set up however. This is how do acceptance testing.

## I understand that I need to add more discussion to this section. There are several details which are useful as time goes on, but not to "get up and running". 12/2002.

I am working with Mark Rivers of BIOCARS/14. He and I discussed how to get the data into his EPICS system. Thus we hope you can just obtain his EPICS drivers. Mark is

using EPICS to control the feedback of a piezoelectric crystal in the monochromater, for beam stability. Talk to him about his EPICS PID controller record.

For the past 3 years, Steve Wasserman does all this through SPEC.

We now discuss details of software setup, emphasis on the register mode.

#### 5.1 SOFTWARE REQUIREMENTS FOR REGISTER MODE.

**5.1.1 DATA OUT TO THE ADCMOD2 MODULES FOR SETUP.** The user must supply initializing information to the ADCMOD2 electrometer. For register based readout, this data sets the electrometer's feed back capacitor gain range, the period for the charge collection measurement, some details of for the DDC112 integrated circuit's operation related to continuous or non-continuous modes, and a go-signal to begin to take data. (For dc storage scope mode, this setup data must also indicate how many samples are required before terminating scan.)

This initializing need only be done upon power up. If you are happy with the gains then it never needs to be repeated. This is a common occurrence with the 20-bit ADC with  $a10^6$  dynamic range. All ADCMOD2 modules are set up identically (assuming they are cabled in, and powered up.)

Internal to the system, over the fiber, between FIBER2 and ADCMOD2, data is passed as serial bit streams 48 bits long. A command to the ADCMOD2 is 48 serial bits long. For the VME host to command the ADCMOD2 it must write three 16 bit hex words to three specified "mailbox" addresses, then write a fourth word to a fourth address to trigger the serial burst. This process is then looped to send additional 48 bit serial bit commands. The specific addresses here can be set in the firmware, and need to be discussed. This system was implemented to give some flexibility – part of the 48 bit command is a device address. Most users at present will only have one TYPE of device out on the fiber link – a varied number of ADCMOD2's. However sector 32 has additional devices.

The following is a basic script. Only 5 commands are really needed: range, pulse, period, conv, go. There is no real order, but typically the go-command comes last. As time goes on I will try to cover some of the "bells and whistle" commands, they are often for de-bugging and test.

Read the following script as AAAADDDD #comment.

For example **f018a000 2range** from the table below means to my system:

Write A24,D16 by PIO transfer to the VME bus. Write to address 0x00f018. The data is 0xa000. The upper data lines (VME supports 32) are don't-care. The upper address lines (VME supports 32 bit addressing) are don't care. VME supports address modifiers

#### Electrometer Users Guide 12-17-02.doc

(AM5:0), and various signals to say the data is a LONG, which style ENDIAN in use etc. While I bring these signals into the Altera chip, as of now I do not decode them. Likely this will be added, as I understand any issues users have with VME bus contention – it is not difficult.

Commands below with comments 2range, 3range, 4range combine to form the 48 bit stream internally. The 5send command is simply are VME write to address 0x00f020 and triggers the internal serial stream.

Thus in the example the range command is (hex) a000 0001 0000

the pulse command is a000 0006 100f

the period command is a000 0007 ffff

the conv command is a000 0005 0280

the go command is a000 0004 0001.

This test script is an input to my C++ code. Yours may vary slightly. It is nice to have the comments out to the side.

40 testdc2 f018a000 2range f0100001 3range f0080000 4range f020cccc 5send f018a000 6pulse f0100006 7pulse f008100f 8pulse f020cccc 9send f018a000 10period f0100007 llperiod f008ffff 12period f020cccc 13send f018a000 14conv f0100005 15conv f0080280 16conv f020cccc 17send f018a000 18go f0100004 1990 f0080001 20qo f020cccc 21send

#### 5.1.2 DATA COMING BACK FROM THE MODULES:

Internal to the system, and transparent to the user, data returning from and ADCMOD2 is again packaged into 48 bits. These bits contain diagnostic, address and parity information and all the data bits from the ddc112. There are dip-switches on each ADCMOD2 printed circuit board. The switch pattern is also embedded in these 48 bits.

Under normal operation only the ddc112 data is presented into the VME accessible registers. There are diagnostic modes. (need to discuss someday...)

We now discuss this VME returned data.

The DDC112 output values go from0x00000 to 0xfffff. (see its data sheet) This entire 20-bit value is given to PASSTH4. But since the VME data word is only 16 bits, we must somewhere loose 4 bits. If there is no signal averaging we simply drop the 4 least significant bits. I note that the bottom 4 bits (16 counts) are noise anyway at this point (16 counts correspond to about 16 pA, the present noise level is above this.) Hence the VME registers see the most significant 16 bits from the 20 bit amplifiers.

At the very beginning I'll keep signal averaging off. When it comes on, the following will apply: There must be a match between which bits are given to the user and the number of signals averaged, otherwise there will be a lot of wasted bits. For example if 8 values are added, then in the sum the least significant 3 bits will have little meaning. I need to work on this but the system will eventually give the most significant bits, and just exactly which bit is MSB will depend on the averaging. Averaging is done in a 30-bit adder so we can add about  $2^{14}$  counts, i.e. 4096 counts.

**The DDC112 has an internal DC offset.** When no current is going into the module, the output value is 0x01000. I tend to see this a little low, e.g. 0x00f7a or thereabouts. You can get a quick indication of the system noise by watching the variations here.

To convert digital counts to actual current, you must know the gains. Typically I do this the easy way and inject a known current from a Keithley source. Again 1pA is about one analog-digital unit (ADU) when range = 000.

## 5.2 SOFTWARE REQUIREMENTS FOR STORAGE SCOPE MODE.

Discussed in the appendix for now.

## 5.3 DEBUGGING MODES – HOW TO READ DIAGNOSTIC BITS INTO THE VME HOST COMPUTER.

I'll expand on this as time goes by. But the idea is that we can make changes to the system to quickly access its health. A typical thing to do would be to replace electrometer data with bits showing which ADCMOD2 is which, that the parity is ok, pick out the dip-switch information etc.

## (6) INSTALLING A SYSTEM AT AN APS BEAMLINE.

**6.1 COSTS** As I understand my job per APS management, our primary mission is to support the CAT beamline personnel. This electronics is thus made available for minimal costs of manufacture about \$2000 - \$3000 total, for the 4- to 16-channel system. This breaks down to \$1000 per VME board, and about \$500 each for an ADCMOD2 board. Thus the incremental cost of adding four more electrometers is \$500. This price is expected to drop more as quantities rise.

**6.2 HARDWARE MAINTENANCE/ SUPPORT**. We in APS will keep your system working on into the future. This is not a trivial undertaking. We have a small staff, but it is more than one person. Our hope is that a number of beamlines and sectors will standardize on how to read the data out, and this will help a lot with configuration management, and keeping spare parts. Support will **not** be 24 hours/day, 7 days per week, but more like "working hours". The best way to have a reliable system is to have some backup channels. It will take us a while to get there.

We are currently installing a system on about 11 beamlines, so there is some comfort in numbers.

We follow ASD quality control and documentation guidelines, again not trivial.

#### 6.3 INSTALLATION ISSUES.

#### FIBER OPTICS.

Each adc module communicates over a pair of fibers. The fiber opticcable comes in fibers to a cable "4-fiber cable". This cable isbreakout fiber, orange color. I supply it.

- (1) How long is your length of fiber?
- (2) How many 4-fiber cables do you wish pulled now?

(3) will your technicians pull it, or do you want ASD techs to help (at \$70 an hour I think)?

(4) Do you wish ASD to terminate your fibers with the required ST-2 connectors?

#### PIN DIODE HOLDERS.

I need to work with each sector to mount the PIN diodes, and make sure the wires are brought out. Then if it is different I can fabricate a new PCB (cost about \$400 qty 10). Please get me this mechanical design info, or use the existing ones.

I already have several designs for the holder for the PIN diode: See sector 2 for an Autocad drawing for two designs (one worked out with John Quintana sector 5 which precisely fits inside a KF-50 cross. I call it the diode6 design.)

Also you can notice a sketch of a round pcb 1.54" diameter. I call it the diode5 design.

I have another fairly complicated design which slides into a KF50 cross like a cylinder in a pipe. Call it diode4.

I will work with people who need ultra high vacuum. Typically they will guide me on how this can be implemented. Otherwise people hold the diodes on with superglue, 5 minute epoxy etc.

(5) Do you wish to work with me to hold the PIN diodes, if so how --standard product or custom?

#### FOIL HOLDER.

Using diode5 or diode6 design, some sectors just want to have another pcb mounted to the irst, and onto it glue down the foil. These "blanks" have no diodes on them, but just simply hold the foil. (FIG DIODE7.dwg) They can be mounted on stand-offsontodiode5 or diode6. Typical stand off distance is 1 cm.

(6) Do you wish to work with me to hold the foils? If so how?

(7) Will you get your own foils? I can tell vendors, but typically I don't worry about the actual foil. I learn much of this from Randy Alkire at sector 19. Other sectors will coat kapton with metal.

#### WIRING TO PIN DIODES.

I can supply the hookup wires from the pindiodes to the adc module, though it is just ordinary wire. I have a minor preference for wire colors (black, brown, red, orange) to try to "standardize". At sector 32 I was able to tightly twist 2 wires together and get good noise levels even with up to 2 meters of twisted pair, something I did not expect at the time. Of course you can also split out, and run 2 pin diodes at one place, 2 at another etc. The PIN diodes are completely floating. Connectors are always a concern. For the designs consider "standard", I have a stock of connector. Don't be surprised if we have to order stuff from Newark.

#### PIN DIODES AND LIGHT.

The pin diodes respond very well to visible light. For example if you have a vacuum gauge, a Bayart -Alpert vac-ion gauge, it will have an orange glowing filament. If the pin diode can see it, it will happily put out lots of false signal.

POWER SUPPLY.

The adc modules require power. Typically this is about 10 VDC, about .4Amps per adc module. The exact voltage is not critical, as I clean it up. People might want to start out with a lab supply such as a BK Precision 1760 \$600 from Digikey, which can run 3 adc modules.

I supply the connector which plugs into the adc module, to bring power into the adc module. I bring this connector out to 2 wires (hot/gnd). From there it is your problem.

Note: as 2003 progresses, and I see that the systems are working well for you all, I will likely build a ping-pong battery system. This system would be charging one small lead acid 12 vdc battery while draining another. Batteries are great for noise immunity. I cannot find a commercial device which will do this. I have a prototype of this, it is about 12" x 8" by 6" or so. One can get lead acid 12 v batteries from newark, the ones I use typically will run a module for about 12 hours.

(8) Do you have your own dc power, is this all ok?

(9) Would you be interested in a battery system, cost approx. \$400, next year ?

#### TESTING.

I'll bring a cart out with a PC, running windows 98, and my own C++ code. I'll put the cards my VME crate, along with a LED bus monitor card. (FIG 6) I'll show you a working system. Then we have to work together to get into your software. I display the data in IDL, and use IDL to compute noise etc.



Fig 6. Typical test stand consists of a VME rack and a PC running Windows. The software for my own testing is simply MSVC++ and IDL (which are linked to pass data from one to the other.)

## 7.0 ISSUES RELATED TO OPERATING THE SYSTEM

Taken from instr4.txt: version 9/29/00, work in progress. Especially need a better explanation of the LED's.

7.1 REBOOTING THE SYSTEM -- TURN ON SEQUENCE.

(1) Turn off the vme crate, the adc modules and any other modules on the fiber chain (e.g. encoder communications module) (nNte you do NOT have to turn off the Heidenhain encoder module)

(2) Turn ON the VME crate.

(3) Reset the PASSTH4 card by pushing the button on it.

(4) Turn on the encoder communcications module. 3 led's come on.

(5) Turn on the ADC modules. 3 led's come on.

(6) Push the Passth4 "reset" button

#### 7.2 DETAILED EXPLANATION OF THE COMMANDS TO THE ADCMOD2

Set up all of the pindiode ADC modules via the fiber optic communication line, to provide the correct data, and to then command the ADC modules to begin taking data.

I write informative comments after a ! symbol. Delete them for your actual script. All numbers are hexidecimal.

#### 7.2.1 RANGE

a000	2range	!command #1 to adc module	
0001	3range	!command 0001 sets the adc's gain rar	ıge
0000	4range	gain being set to 0000	

Allowed values for adc gain range:

a000 0001 Orrr set the gain of the ADC unit. rrr can be

000	least gain	range , refer t	to it as gain = 1
001	gain = 220	/(1*12.5) = 17	7.6
010	gain = 220	/(2*12.5) = 8.	.80
011	gain = 220	/(3*12.5) = 5.	.87
100	gain = 220	/ (4*12.5) = 4.	. 40
101	gain = 220	/(5*12.5) = 3.	.52
110	gain = 220	/(6*12.5) = 2.	.93
111	gain = 220	/(7*12.5) = 2.	.51

With gain range set to 000 (command a000 0001 0000), and with the conversion time set to 0200 (command a000 0005 0200) and with 10 nA going into the system, the counts will read 14260 (approx). There is a baseline of approximately 4096 counts. So 10 nA creates a signal of 14260 - 4096 = 10164 counts in 820 us. So 1 pA would create just about 1 count.

7.2.2 NUMBER OF DATA POINT	S BACK FROM ADCMOD - digital storage scope
mode only	
a000 6pulse	!command #2 to adc module
0006 7pulse	!command 0006 sets number of data points
1000 8pulse	here we ask for 1000 data points!
7.2.3 PERIOD	
a000 10period	!command #3 to adc module
0007 11period	!command 0007 is currently not supported
ffff 12period	so leave this exactly as shown!

7.2.4 CONVERSION TIME - conv command a000 14conv !command #4 to adc module 0005 15conv !command 0005 sets conversion time 0280 16conv !here we set to 0280 (comment 5) So handy rule of thumb: With conversion time 0200 hex (820us), range = 000, 1 count = 1 pA. Command a000 0005 xxxx sets the conversion time (for a "ping", "pong" is the same). This number, specified by the xxxx is converted to seconds as follows: microseconds = xxxx (hex --> decimal) \* 1.6 + 0.6 example: xxxx = 0280 0280 in hex =  $2*16^2 + 8*16^1 + 0$  $= 2 \times 256 + 8 \times 16 + 0$ = 620 decimal time (microseconds) = (620 \* 1.6) + 0.6= 992.6 The allowed range for this parameter is in hex 2000 >= xxxx >= 180 (hex) thus the conversion time is presently limited to 13.107 milliseconds >= conversion time >= 615 microseconds. Thus the fastest data currently comes at 1/615 microseconds = 1.6 kilosamples/second.

#### 8.0 DESIGN FILES.

This section will grow into a description of the design files - OrCAD, printed circuit board, Altera and AutoCAD. For now know that there is a block diagram on the ESG shared disk. We have bills of materials into excel spread sheets. The altera designs have been extensively simulated.

#### 9.0 FUTURE FEATURES.

Much will depend on feedback from users. This is just a laundry list.

Rechargeable ping-pong battery supply (while one is draining, the other is charging, the ADCMOD2's never see line current).

Better diagnostics.

24 of 38

Make sure the LED's are useful to the user.

### 10.0 APPENDIX. DEBUGGING INFORMATION. USER DOES NOT NEED THIS.

EMBEDDED DE-BUGGING DATA. NOT DISPLAYED IN REGISTER MODE. 1 adc module, 1 encoder communciation module a-a-d- a-a-d- a-a-d a-a-d ping pong i.e. quantity 24 32 bit words going into memory per scripted pulse See also "Interpretation of the Data, comment 6. The 8 a- 's listed just now each are tagged by 8 unique tags. All would come in from the same fiber optic channel. 2 adc modules, 1 encoder communciation module a-a-a-d- a-a-a-d- a-a-a-d- a-a-a-dping pong i.e. quantity 40 32 bit words going into memory per scripted "pulse" 3 adc modules, 1 encoder communciation module ping pong i.e. quantity 56 32 bit words going into memory per scripted "pulse" \_\_\_\_\_\_ C. EXAMPLE SCRIPTS: THESE SCRIPTS ARE WRITTEN IN THE FORMAT FOR SECTOR 32. THEY ARE OF LIMITED VALUE TO OTHER SECTORS, BUT HELP ME DEBUG. Note #1: In these scripts, the number e.g. 0100 on the first line, is to be sent to vme location b002 Note #2: The English comment just afterwards is a human oriented comment, not to be sent out. Note #3: The line 0000 ----GO-GO really means send anything (0000 is ok) to vme location b004 to execute the commands loaded into b002. Example 1: Goal: obtain just a little bit of data, to quickly read the encoder etc. 0100 ----setupacq 0000 commandmode 0000 val 0000 val 0000 val 0000 resv 0001 heartbeat 0000 timeout 0100 blksize32 !very short block 0001 numblks 2000 VSBaddr2000 0020 VSBaddr0020

0000 blkincr

0100 blkincr 2000 VSBstatus 0000 status--0200 ----remotecommand 0003 commandmode 0000 zero 001a size+4 0000 lleading zero a000 2range 0001 3range 0000 4range 0000 5break a000 6pulse 0006 7pulse 0100 8pulse 0000 9break a000 10period 0007 11period ffff 12period 0000 13break a000 14conv 0005 15conv 0280 16conv 0000 17break a000 18flush 0001 19flush 0007 20flush 0000 21trailzero 0000 22trailzero 0300 ----GiveData 0000 val 0000 val 0000 val--0000 ----GO-GO \_\_\_\_\_ Example 2: Goal: once you have been taking data for awhile, and you ONLY want to change the gain range of the ADC: 0200 ----remotecommand 0003 commandmode 0000 zero 0009 size+4 !note size + 4 changed. 9 lines in this command 0000 1leading\_zero a000 2range 0001 3range 0003 4range !new range, now 011 (3 in binary) 0000 5break 0300 ----GiveData 0000 val 0000 val 0000 val--

0000 -----GO-GO !sent to vme address b004 -----Example 3: Goal: once you have been taking data for awhile, and you ONLY want to change the integration time of the ADC: 0200 ----remotecommand 0003 commandmode 0000 zero 0009 size+4 0000 1leading\_zero a000 2conv 0005 3conv 0200 4conv !new conversion time, hex 0200, or 820 us 0000 5break 0300 ----GiveData 0000 val 0000 val 0000 val--0000 -----GO-GO !sent to vme address b004 \_\_\_\_\_ Example 4: Goal: Using the adc module in non-continuous mode. I have seen that this will allow the unit to process up to 5.8 uA of current without saturation.Conversion input of hex 24 yields a conversion time of1.6 \* (0x24 --> decimal 36) --> 58.2 microseconds. Period input of 300 yields a period time of 1228.6 microseconds (per L092900.D01). (i.e. 300 hex --> 768, 768 \* 1.6 = 1228, 1228 plus 0.6 = 1228.6 microseconds). [ref notes 09/21/2000] 48 dcltest1 b0020100 ----setupacq b0020000 commandmode b0020000 val b0020000 val b0020000 val b0020000 resv B0020001 heartbeat B0020000 timeout B0021000 blksize32c000 B0020001 numblks12 B0022000 VSBaddr2000 b0020004 VSBaddr0020 b0020000 blkincr b0020100 blkincr b0022000 VSBstatus

b0040000 --GO--

b0020200 ----remotecommand

b0020003 commandmode b0020000 zero b002001a size+4 b0020000 1leading zero b002a000 2range b0020001 3range b0020000 4range b0020000 5break b002a000 6pulse b0020006 7pulse b0020800 8pulsefff0 b0020000 9break b002a000 10period b0020007 11period b0020300 12period b0020000 13break b002a000 14conv b0020005 15conv b0020024 16conv b0020000 17break b002a000 18flush b0020001 19flush b0020007 20flush b0020000 21trailzero b0020000 22trailzero b0020300 ----GiveData b0020000 val b0020000 val b0020000 val--B0040000 ----GO-GO

\_\_\_\_\_

D. INTERPRETING THE DATA. LIMITED USE FOR PEOPLE RUNNING REGISTER MODE, MAINLY FOR DIGITAL STORAGE SCOPE MODE. HELPFUL FOR ME TO DEBUG.

I include here a sample readout

of the PINDIODE system, with 2 ADC modules connected. The ADC modules are connected to fiber optic channels 1 and 4. The encoder communcications module is connected to fiber optics channel 3. We can label the columns of data (printed out in c by a 10%x format) as follows:

LONG	ADDR	CONTENTS	COMMENTS
WORD	VME		
NUMBER			

First Note:recall in the tutorial script the lines2000 VSBstatus!address where DC1 status register goes0000 status--!this status register upon success reads 20001which stated that the vme status register should be address 2000 0000

0	20000000	20001	!VME transfer status register
1	20000004	fffffff	! junk
2	20000008	ffffdff	!this area of memory is never

3 4 5 6 7	2000000c 20000010 20000014 20000018 2000001c	fdffffff 7ff7efff fffffffe dffffbef fffffbf	!written to. !end of junk
Second Note 2000 VSBadd 0020 VSBadd going.	e: recall in Ar2000 Ar0020	the tutoria !higher orde !lower order	el script the lines: er address of where data is going e VME/VSB address of where data is
8	20000020	a1ee0c00	lbegin data #1 (first pulse)
9	20000020	11	Inote address 2000 0020
a	20000021	alef0c2f	!was set by these lines
b	2000002c	fffc0014	
C	20000030	alee1c00	
d	20000034	21	
е	20000038	alef1c40	
f	200003c	698a0024	
10	20000040	d0fe0017	
11	20000044	60040013	
12	20000048	alee0768	
13	2000004c	2a2c0031	
14	20000050	alef07ac	
15	20000054	5d880034	
16	20000058	alee1400	
17	2000005c	41	
18	20000060	aleil/cl	
19	20000064	learuu44	
1a 1b	20000008 2000006c	60040023	
10	200000000	00010025	blank line inserted for
clarity			
1c	20000070	a1ee0980	
1d	20000074	29800051	!adc still "settling down"
le	20000078	a1ef0940	-
1f	2000007c	330c0054	
20	20000080	alee1800	
21	2000084	61	
22	20000088	alef18c0	
23	2000008c	7400064	
24	20000090	d0fc0017	
25	20000094	60050033	
26	20000098	alee0000	
27	20000090	12000071	
28	200000a0	aler0180	
29 20	200000a4 200000a9	∠34000/4 a1ee1000	
∠a 2h	200000000	81	
20	200000000	alef1020	
2d	200000b4	5980084	
2e	200000b8	d0fc0017	!encoder is at 176005
2f	200000bc	60050043	!end of pulse 1's data
			!blank line inserted for clarity
30	200000c0	alee0dc0	pulse #2
31	200000c4	3ae00091	
32	200000c8	alef0c80	

	33	200000cc	6£940094	
	34	200000d0	alee1c80	
	35	200000d4	ecd00a1	
	36	200000d8	alef1c40	
	37	200000dc	fec00a4	
	38	200000e0	d0fc0017	
	39	200000e4	60040053	
	3a	200000e8	a1ee0400	
	3b	200000ec	b1	
	3c	200000f0	a1ef0560	
	3d	200000f4	52fc00b4	
	3e	200000f8	alee14a0	
	3f	200000fc	ef000c1	
	40	20000100	a1ef1540	
	41	20000104	10d300c4	
	42	20000108	d0fe0017	
	43	2000010c	60040063	
				!blank line inserted for
clarity				
-	44	20000110	a1ee0980	
	45	20000114	378c00d1	
	46	20000118	a1ef0900	
	47	2000011c	3e0c00d4	
	48	20000120	alee18a0	
	49	20000124	e3300e1	
	4a	20000128	a1ef1820	
	4b	2000012c	fa200e4	
	4c	20000130	d0fe0017	
	4d	20000134	60050073	
	4e	20000138	alee0160	
	4f	2000013c	411c00f1	
	50	20000140	alef01a0	
	51	20000144	3d9000f4	
	52	20000148	alee10c0	
	53	2000014c	db20101	
	54	20000150	a1ef1080	
	55	20000154	ef30104	
	56	20000158	d0fe0017	
	57	2000015c	60040083	
	58	20000160	alee0d80	pulse #3
	59	20000164	3cbc0111	
	5a	20000168	alef0da0	
	5b	2000016c	50100114	
	5c	20000170	alee1c60	
	5d	20000174	fbc0121	
	5e	20000178	alef1c60	
	5f	2000017c	fb20124	
	60	20000180	d0fe0017	
	61	20000184	60050093	
	62	20000188	alee0400	
	63	2000018c	131	!don't worry,I fixed this
later	•			
	64	20000190	alef05e0	
	65	20000194	44480134	
	66	20000198	alee1400	
	67	2000019c	ffc0141	

68	200001=0	a1ef1560
C 0	20000140	10100144
69	200001a4	10180144
6a	200001a8	d0fe0017
6b	200001ac	600500a3
бc	200001b0	alee0940
6d	200001b4	384c0151
60	20000168	alof0920
6 E	200001b0	20700154
01 D	20000100	30/00134
70	20000100	aleel8CU
71	200001c4	e230161
72	200001c8	alef1860
73	200001cc	fb20164
74	200001d0	d0fc0017
75	200001d4	600400b3
76	200001d8	alee0180
77	200001da	42000171
70	20000140	42000171
/0	20000100	aleluluu
/9	200001e4	3e000174
7a	200001e8	alee10c0
7b	200001ec	dd20181
7c	200001f0	alef1020
7d	200001f4	£570184
7e	200001f8	d0fc0017
7f	200001fc	600400c3
9 A	20000210	21000d60
01	20000200	a100000
81	20000204	36300191
82	20000208	aleiud2u
83	2000020c	47480194
84	20000210	aleelc60
85	20000214	fbc01a1
86	20000218	alef1c00
87	2000021c	fc601a4
88	20000220	d0fc0017
89	20000224	600400d3
82	20000221	21000400
oh	20000220	1b1
0	20000220	
8C	20000230	aleruscu
8d	20000234	422c01b4
8e	20000238	alee1540
8f	2000023c	102201c1
90	20000240	alef1540
91	20000244	102301c4
92	20000248	d0fc0017
93	2000024c	600400e3
94	20000210	210010000
	20000250	201-01-01
95	20000254	391C0101
96	20000258	alei09c0
9.7	2000025c	3c5401d4
98	20000260	alee1880
99	20000264	e6301e1
9a	20000268	alef1840
9b	2000026c	f2901e4
90	20000270	d0fc0017
94	20000274	600500f2
90	20000274	21000120
9e		A1400151
ЭI	2000027C	414UULTI
aU	20000280	alet01a0

al	20000284	3d5c01f4
a2	20000288	alee1080
a3	2000028c	df10201
a4	20000290	alef1020
a5	20000294	£620204
аб	20000298	d0fe0017
a7	2000029c	60050103
a8	200002a0	alee0d40
a9	200002a4	3e500211
aa	200002a8	alef0de0
ab	200002ac	44500214
ac	200002b0	aleeld40
ad	200002b4	10100221
ae	200002b8	alef1c60
af	200002bc	£470224
b0	200002c0	d0fe0017
b1	200002c4	60040113
b2	200002c8	alee0400
b3	200002cc	231
b4	200002d0	a1ef0560
b5	200002d4	3e840234
b6	200002d8	alee1540
b7	200002dc	104d0241
b8	200002e0	alef1460
b9	200002e4	£770244
ba	200002e8	d0fe0017
bb	200002ec	60040123
bc	200002£0	alee0920
bd	200002£4	392c0251
be	200002f8	alef0900
bi	200002fc	38c00254
CU	20000300	alee18c0
CI	20000304	e430261
c2	20000308	alei1880
c3	2000030c	e550264
c4	20000310	dufeuul7
C5	20000314	60040133
C6	20000318	alee0160
C /	20000310	41000271
C8	20000320	aler0160
09	20000324	38440274
ca	20000328	aleeludu aafoool
CD	20000320	uuruz81
CC	∠0000330	aleriusu
ca	∠∪∪∪∪334	eucuz84

etc. etc. etc

#### HOW TO INTERPRET THE DIGITAL STORAGE SCOPE MODE DATA.

I have not finished de-bugging the parity, so don't worry about it yet.

Take the last adc data given above:

alef 1080 e0c 0284

```
parse it:
a = comes from the binary pattern 1010, indicates an adc module
1 = comes from the binary pattern 0001, ignore this
e comes from a binary pattern 1110
       in general nnnn
       nnnn = 1 1 1 test, test usually is 0, so nnnn = 1110
       the test bit goes hi when we run a specific diagnostic on the
adcmod,
       the "Burr Brown DDC112 ADC test" diagnostic.
       See DDC112 data sheet Figure 8.
        (When test is turned on by a000 000f vvvv, when vvvv = 7,
       adc reads 300,000 decimal. when vvvv = 3, adc reads 150,000
decimal.)
       clear test with a000 0002 0000.
f comes from a binary pattern 1111
       in general pppp
       pppp = switch switch switch
            = switch settings on the adcmod (a tag for this particular
            module, set by the human using switches)
1 comes from a binary pattern 0001
       in general qqqq
       qqqq = range range tag
       e.q. if you were using range 000 then this should be ggg = 000
       e.g. if you were using range 111 then this should be qqq = 111
       see Burr Brown DDC112 data sheet table I, page 7.
       the last bit, the 4th q bit "tag", indicates which of the two
DDC112
       burr brown chips on the printed circuit board is being read
       0=DOUTA=U5, 1=DOUTB=U8 (see comment 6)
0 comes from a binary pattern 0000
       in general ssss
       ssss = 1st/2nd_readout ping/pong parity19:16 parity15:12
       (see comment 6 below for explanation of 1st/2nd, ping/pong,
       as these are the second and third "tag" bits)
8 comes from a binary pattern 1000
       in general tttt
       tttt = parity11:8 parity7:4 parity3:0 notused
Actual data comes next and is 00e0c note the missing leading 0 in
my printout, just an artifact of my c code. Note that the baseline
of the adc converter is 01000 (see Burr Brown data sheet, Table VII
for more details of data format). Hex 01000 is decimal 4096.
next comes
28, the (hex) 28th data coming in
next comes
4 data coming in from fiber optic channel 4
_____
comment 6
```

#### Electrometer Users Guide 12-17-02.doc

The tags for which adc readout is which. If you put together the 4th "q" bit, with the 1st and 2nd "s" bit you get 3 bits. These tag bits come out in the following order: tag bits: 011 111 001 101 010 110 000 100 The interpretation of these 8 permutations of 3 bits are: DOUT is pin 16 of either the A or B adc integrated circuit u5,u8 Input 2 is pin 28 of the Burr Brown DDC112, Input 1 is pin 1. J1 is the connector at the edge of the printed circuit card. For explanation of conv look at Figure8 of DDC112 data sheet. tag altera DDC112 J1 wire color bits pin pin pin 011 DOUTA input2 conv=hi u5.28 IN2A J1.3 RED ping PINDIODE2 111 DOUTB input2 conv=hi u8.28 IN2B J1.7 YELLOW PINDIODE4 ping 001 DOUTA input1 conv=hi u5.1 IN1A J1.1 BROWN PINDIODE1 ping 101 DOUTB input1 conv=hi u8.1 IN1B J1.5 RED PINDIODE3 ping 011 DOUTA input2 conv=lo u5.28 IN2A J1.3 RED PINDIODE2 pong 111 DOUTB input2 conv=lo u8.28IN2B J1.7YELLOW001 DOUTA input1 conv=lo u5.1IN1A J1.1BROWN PINDIODE4 pong PINDIODE1 pong 101 DOUTB input1 conv=lo u8.1 IN1B J1.5 RED PINDIODE3 pong

So when you for example pull out the tag pattern 011 from the "q" and "s" bits, you have read the "pong" read of PINDIODE2, coming to you over a red wire.

Note that when we look at the printed circuit board containing the pindiodes, looking at the diode side of the board, we see them in the order

#### PINDIODE1

PINDIODE2 hole PINDIODE4

#### PINDIODE3

x-rays are coming out of the page. Pindiode 1 is physically on top.

PARITY. A modulo two sum of bits, useful to detect bit errors. Not yet fully de-bugged.

For example parity19:16
count number of 1's in the bit pattern.
If even, parity bit = 0

IF odd, parity bit = 1 example1, (hex) data is 
 o
 8
 8
 hex

 1000
 1000
 1000
 1000
 binary

 19:16
 15:12
 11:8
 7:4
 3:0
 position

 1
 1
 1
 1
 number of 1's

 1
 1
 1
 1
 parity
 88888, that is example2, (hex) data is cccc, that is ccccchex11001100110011001100binary19:1615:1211:87:43:0position22222number of 1's00000parity example3, (hex) data is 5710f, that is 5 7 1 0 f hex binary position 0101 0111 0001 0000 1111 19:16 15:12 11:8 7:4 3:0 2 3 1 0 4 0 1 1 0 0 number of 1's parity 0 \_\_\_\_\_ E. HARDWARE SYSTEM. E.1 block diagram. <--> DC1 <--> PASSTH <--> FIBEROPTIC <--> ADC SBS module Bus front Master vme VSB <--> Encoder backplane backplane grey Communication cable Module <--> Chrislin VSB Memory VME backplane Card Set up commands come from the user's computer via the SBS VME master, and are passed in a standard manner, over the VME backplane to the commercially supplied DC1 card (from Access Dynamics 505 345 7637). This is accomplished by a VME A16/A24 D16 read/write to the DC1's "mailbox" address, hex b002. The scripts listed above are loaded, one word (16 bits) after another, into the mailbox address. When the user wishes the DC1 to take action based on this script, the user simply writes any word into address b004, which triggers an interrupt on the DC1,

and forces it to execute the script.

The DC1 interprets what it receives. This is covered in some detail in

the manual "Access Dynamics Inc. Electrometer Interface and Firmware". In essence there are 3 types of commands, which are decoded from the scripts: o DC1 Static Setup Commands (we called these Task 1 commands above), o DC1 Dynamic commands (not covered yet in this manual, future...), and 0 DC1 remote commands (we called these Task 2 commands above). The DC1 eventually must pass data to the Chrislin memory. Hence this memory must be setup to have a proper VME/VSB address which matches that in the DC1 script. E.2 MASTER ADC MODULE. The ADC modules are set up by Task 2-type commands to the DC1, which in turn passes a serial bit stream to PASSTH and on to FIBEROPTIC card, and on to the (multiple) ADC modules possibly connected in parallel with the encoder communciation module. At the present time all Task 2 remote commands are addressed to the adc module. Thus we must discuss how the encoder module is commanded/queried. One ADC module is the "master" module. When data is returned from the master module to the VME fiber optic module, a request is sent (by PASSTH) out to all modules requesting Encoder Communciation Module data. No master module, NO requests to the Encoder. There need not be an Encoder Communications module in the system. In reality every other data returned from the master generates an encoder request. Hence the returned data will have patterns such as a-a-d, where the a-'s are returned data from the master. THE MODULE CONNECTED TO FIBER CHANNEL 4 IS THE MASTER MODULE. An infinite loop would be created if the Encoder Communciation Module were connected here! E.3 the LED's. LED's on the DC1 card Underneath the label "serial port" are 4 small green LED's DC1 LED1 : DC1 wants data DC1 LED2 : DC1 is getting and processing data DC1 LED3 : DC1 LED4 : DC1 sending out serial data (remote command being executed. The DC1 requires the amount of data specified in the block size (times the number of blocks if this is greater than 1). As the DC1 receives 32 bit words, it decrements the quantity. When this reaches 0, the DC1 is

satified, and turns off LED1 and LED2. If LED2 goes off to early, LED1

will remain ON, indicating the DC1 requres more data. The DC1 however will flush whatever data it has received so far on to the Chrislin memory.

LED's on PASSTH VME card.

PASSTH LED1 indicates that the crystal clock is running. This LED should come on upon power up. It is nearest the VME card horizontally. PASSTH LED2 indicates that the PASSTH card is receiving data from the FIBER OPTIC card. Data is coming into the system from SOMEWHERE, either an adc module, or the encoder communciations module. It is mounted to the right of the first LED.

LED's on the FIBEROPTIC vme card. Front LED, next to PCB: Good to Go Front LED, away from PCB: Receiving A's data Back LED, next to PCB: Handshake of data passing from Fiber2.pcb to passthl.pcb Back LED, away from PCB: Fiber2 putting out data to adcmod's etc

LED's on the ADC module or Encoder Communcication Module. 4 Red LED's. ???? I may have the order wrong...

LED ADC3 LED ADC1 LED ADC4 LED ADC2 LED ADC1 = ON when there is a bit stream coming into the module It is a stated fact about fiber optic communications that this bit stream may or may not be carrying any information. However it must ΒE THERE, or there is no chance of information flow. LED ADC2 = ON when there is a bit stream exiting the module. It is a stated fact about fiber optic communications that this bit stream may or may not be carrying any information. However it must BE THERE, or there is no chance of information flow. LED ADC3 = ON when there is information (1's actually) coming into the ADC module from the fiber optic module. LED ADC4 = ON when the ADC module Altera integrated circuit is delivering a clock to its fiber optic transmitter integrated

Thus upon power up, LED's ADC1,2,4 should be ON. This assumes that the fiber cables are connected, and it assumes that the fiber optic module is working back at the vme rack.

circuit.

DIP switches on fiber2.pcb switch 1 up indicates that this fiber2 board is the master fiber2 board. Data coming in on channel 0 of the master fiber2 pcb will result in the

b002a000 2range

passthl board sending out commands to read the encoder. switch 2 up indicates that fiber2 is putting out test data, a pattern of a000 0000 0000's to the modules. This is useful for debugging but should be turned off (switch down) during normal operation switch 5 up to enable data from channel 4

switch 6 up to enable data from channel 3 switch 7 up to enable data from channel 2 switch 8 up to enable data from channel 1 note that channel 0 is always enabled. DIP switches on ADCMOD sw 1 always leave up (on) sw 2 when "off" puts out test pattern of a000 0000 0000 sw 3 lsb of tag sw 4 tag sw 5 tag sw 6 msb of tag e.g. up up up up dn dn would be normal operation of a module, with the tag or module address being "c". We would get data back looking like "alec xxxx xxxx" etc. APPENDIX 1. This is from the DCl script I used for many tests dcl\_sv2.dat. 48 numberlines b0020100 ----setupacq b0020000 commandmode b0020000 val b0020000 val b0020000 val b0020000 resv B0020001 heartbeat B0020000 timeout B0028000 blksize32c000 B0020001 numblks12 B0022000 VSBaddr2000 b0020020 VSBaddr0020 b0020000 blkincr b0020100 blkincr b0022000 VSBstatus b0020000 status-b0040000 --GO-b0020200 ----remotecommand b0020003 commandmode b0020000 zero b002001a size+4 b0020000 1leading\_zero

b0020001 3range b0020000 4range b0020000 5break b002a000 6pulse b0020006 7pulse b0021000 8pulsefff0 b0020000 9break b002a000 10period b0020007 11period b002ffff 12period b0020000 13break b002a000 14conv b0020005 15conv b0020280 16conv b0020000 17break b002a000 18flush b0020000 19flush b0020000 20flush b0020000 21trailzero b0020000 22trailzero b0020300 ----GiveData b0020000 val b0020000 val b0020000 val--B0040000 ----GO-GO